151712 160127

US

For US and Canadian markets only.
For other markets see last page.

HemoCue® Hb 201 Microcuvettes

The HemoCue Hb 201 Microcuvettes are designed for use with the HemoCue Hb 201⁺ Analyzer and the HemoCue Hb 201 DM Analyzer (referred to as the HemoCue Hb 201 Analyzer in this document). HemoCue Hb 201 Microcuvettes are available in individual packages or in vials. Please read the relevant Operating Manual for proper use of the system¹.

Intended Purpose/Intended Use

Quantitative determination of hemoglobin in capillary, venous and arterial whole blood, using a specially designed analyzer, the HemoCue Hb 201 Analyzer, and specially designed microcuvettes, the HemoCue Hb 201 Microcuvettes. HemoCue Hb 201 Microcuvettes are for *In Vitro* Diagnostic use only. The HemoCue Hb 201 Analyzer is only to be used with HemoCue Hb 201 Microcuvettes.

IVD Medical Device Directive

The HemoCue Hb 201 Microcuvettes comply with the IVD Medical Device Directive 98/79/EC and carry the CE mark

Principles of the method/procedure

Principle of the method

The reaction in the microcuvette is a modified azidemethemoglobin reaction. The erythrocytes are hemolyzed to release the hemoglobin. The hemoglobin is converted to methemoglobin and then combined with azide to form azidemethemoglobin. The measurement takes place in the analyzer in which the transmittance is measured and the absorbance and hemoglobin level is calculated. The absorbance is directly proportional to the hemoglobin concentration.

Principle of the procedure

The system consists of an analyzer together with microcuvettes. The microcuvette serves both as a pipette and as a measuring cuvette and is for single-use only. A blood sample of approximately 10 µL is drawn into the cavity by capillary action. The analyzer measures at two wavelengths in order to compensate for turbidity, and the hemoglobin level is calculated and presented. The HemoCue Hb 201 system is calibrated against the international reference method for hemoglobin determination, ICSH³ and needs no further calibration.

Composition

The microcuvette is made of polystyrene plastic. Reagents; <600 µg/g microcuvette sodium deoxycholate, <300 µg/g microcuvette sodium azide, <300 µg/g microcuvette sodium nitrite, <350 µg/g microcuvette nonreactive ingredients.

Warning and precautions

The microcuvettes are for *In Vitro* Diagnostic use only. Always handle blood specimens with care as they may be infectious. Consult local environmental authorities for proper disposal.

Storage and handling

Use the HemoCue Hb 201 Microcuvettes prior to their expiry date. The expiry date is printed on each package.

Storage for microcuvettes kept in a vial

The microcuvettes are to be stored at room temperature (15–30 °C, 59–86 °F) and in a dry place. Once the seal is broken the microcuvettes are stable for three months. Always keep the vial properly closed.

Storage for individually packaged cuvettes

The microcuvettes are to be stored at room temperature (15–30 °C, 59–86 °F) and in a dry place.

Specimen collection and preparation

Capillary, venous or arterial blood may be used. Appropriate anticoagulants (e.g. EDTA or heparin) may be used, preferably in solid form to avoid dilutional effects. Mix all specimen tubes thoroughly on a mechanical mixer for at least 2 minutes or invert the tube 8-10 times by hand. Hemoglobin remains unchanged for days, provided that the blood does not become infected. If the specimen has been stored in a refrigerator, it will be viscous and should reach room temperature before mixing².


Directions for use

Materials required

- HemoCue Hb 201 Analyzer
- HemoCue Hb 201 Microcuvettes
- Lancet (for capillary samples)
- Pipette or other transfer device (for venous, arterial or control material samples)
- Lint-free wipe (non-fraying)
- Hydrophobic surface (for venous, arterial or control material samples)

Procedure

The operating temperature of the HemoCue Hb 201 system is 15-30 °C (59-86 °F). Please read the Operating Manual for proper use of the system¹. See relevant manual for information on repeat capillary sampling. For further information please contact HemoCue.

1. Optical eye
2. Filling end

Remove the microcuvette from the package. Close the vial. Hold the microcuvette opposite the filling end and bring into contact with the specimen, allowing the cavity to fill completely. Always avoid touching the optical eye. Do not refill the cavity of the microcuvette. Note: The microcuvette should be filled within 3 minutes after the microcuvette has been taken out of its package.

When collection tubes or control material are used, dispense an aliquot of the well-mixed specimen onto a hydrophobic surface. Note: The microcuvette should be filled within 3 minutes after the microcuvette has been taken out of its package.

Wipe off the outside of the microcuvette with a clean lint-free wipe, being careful not to touch the open end. If air bubbles are seen in the optical eye of the microcuvette, discard the microcuvette and take a new sample. Small air bubbles around the edge do not influence the result.

Place the microcuvette into the cuvette holder and start measurement as soon as possible but no later than 10 minutes after filling the microcuvette by gently pushing the cuvette holder to its measuring position. Close the cuvette holder. The measuring time is 15-60 seconds for Hb values below 20 g/dL. Pull the cuvette holder out to its loading position and discard the used microcuvette.

Quality control

The HemoCue Hb 201 Analyzer has an internal electronic self-test. Every time the analyzer is turned on, it will automatically verify the measurement performance. This test is performed at regular intervals if the analyzer remains switched on. Follow local guidelines regarding quality control procedures. If a quality control test is required by local or other regulations and therefore should be performed, only use controls recommended by HemoCue. For further information regarding controls contact HemoCue AB.

Limitations of the method

- The measurement needs to be started no later than 10 minutes after filling the microcuvette.
- Mixing samples for an extended period can produce increased oxygen pressure and viscosity that may give falsely results.
- If “HHH” is displayed, the result exceeds the measuring range of the system.
- Values above 23.5 g/dL (235 g/L, 14.6 mmol/L) must be confirmed using a suitable laboratory method.
- Following substances have not been found to interfere: Acetaminophen (20 mg/dL), ascorbic acid (3 mg/dL), conjugated bilirubin (40 mg/dL), unconjugated bilirubin (20 mg/dL), creatinine (30 mg/dL), ibuprofen (40 mg/dL), leukocytes ($600 \times 10^9/L$), lipemia (intralipid 4000 mg/L, triglycerides approximately 1000 mg/dL), salicylic acid (50 mg/dL), tetracycline (20 mg/dL), trombocytes ($2100 \times 10^9/L$), urea (500 mg/dL), uric acid (20 mg/dL). The highest concentration or percentages tested is referred to in brackets. Interference studies have been performed according to NCCLS Document EP-7⁴.
- pH values between 6.3-9.0 do not interfere with the system.
- Sulfhemoglobin is not measured with this method.

Expected values²

Children	11.0-14.0 g/dL (110-140 g/L, 6.8-8.7 mmol/L)
Women	12.0-15.0 g/dL (120-150 g/L, 7.4-9.3 mmol/L)
Men	13.0-17.0 g/dL (130-170 g/L, 8.1-10.5 mmol/L)

Children, 2 years to teenage, gradually increase to adult values.

Due to a wide range of conditions (dietary, geographical, etc) which affect normal values, it is recommended that each laboratory establish its own normal range.

Specific performance characteristics

Within-run and Total precision

Within-run precision was determined according to the NCCLS Document EP5-A⁵. The results given below in “Within-run precision” and “Total precision” come from 1 batch of HemoCue Hb 201 Microcuvettes and 5 HemoCue Hb 201⁺ Analyzers. No recalibration was performed during the analyzing period. Commercially available controls at 2 different levels were used. The hemoglobin concentration was measured in duplicate twice a day, morning and afternoon, during 20 consecutive days.

Control level	N	\bar{x} g/dL	Within-run Precision		Total Precision	
			SD g/dL	CV %	SD g/dL	CV %
1	400	7.70	0.057	0.74	0.100	1.30
2	400	15.33	0.078	0.51	0.109	0.71

Accuracy

The results of the comparison studies between the HemoCue Hb 201 system and the International Council for Standardization in Haematology method (ICSH)³ are summarized in the table below. The study was performed on four HemoCue Hb 201⁺ Analyzer which had not been recalibrated during the study period.

Study	N	Min g/dL	Max g/dL	Regression line	Correlation coefficient (r)
1	498	4.1	20.5	$Y = 1.009X - 0.008$	0.998
2	103	9.0	17.7	$Y = 0.971X + 0.584$	0.934

1 = ICSH (Cyanmethemoglobin method), venous EDTA blood, multicenter study

2 = ICSH (Cyanmethemoglobin method), capillary blood

Bibliography

1. HemoCue Hb 201 Operating Manuals
2. Dacie and Lewis, Practical Haematology, Ninth Edition, 2001.
3. Reference and Selected Procedures for the Quantitative Determination of Hemoglobin in Blood; Approved Standard NCCLS Document H15-A
4. Interference testing in clinical chemistry NCCLS approved guideline; NCCLS Document EP7
5. Evaluation of precision performance of clinical chemistry devices; approved Guideline NCCLS Document EP5-A

Symbols used

	Caution		Batch code		Use by (Year Month Day)
	In vitro diagnostic medical device		Catalog number		Date of opening
	Do not reuse		Temperature limitation		Open vial expiry date. Must not exceed the "Use by" date
	CE mark		Consult instructions for use		
			This product is covered by one or more patents		

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Bibliographie

1. Manuels d'utilisation HemoCue Hb 201
2. Dacie and Lewis, Practical Haematology, Ninth Edition, 2001
3. Reference and Selected Procedures for the Quantitative Determination of Hemoglobin in Blood; Approved Standard NCCLS Document H15-A
4. Interference testing in clinical chemistry NCCLS approved guideline; NCCLS Document EP7
5. Evaluation of precision performance of clinical chemistry devices; approved Guideline NCCLS Document EP5-A

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1 = sânge venos EDTA, studiu multicentric
 2 = sânge capilar

Bibliografie

1. Manual de utilizare HemoCue Hb 201*
2. Instrucțiuni de utilizare HemoCue Hb 201 DM (Manual de utilizare)
3. Dacie and Lewis, Practical Haematology, Ninth Edition, 2001
4. Reference and Selected Procedures for the Quantitative Determination of Hemoglobin in Blood; Approved Standard; CLSI Document H15-A
5. Interference Testing in Clinical Chemistry; Approved Guideline; CLSI Document EP7
6. Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline; CLSI Document EP5-A

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1. Manual de funcionamiento de HemoCue Hb 201
2. Instrucciones de uso de HemoCue Hb 201 DM (Manual de funcionamiento)
3. Dacie and Lewis, Practical Haematology, Ninth Edition, 2001
4. Reference and Selected Procedures for the Quantitative Determination of Hemoglobin in Blood; Approved Standard; CLSI Document H15-A
5. Interference Testing in Clinical Chemistry; Approved Guideline; CLSI Document EP7
6. Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline; CLSI Document EP5-A

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
Bibliographie

1. Manuel d'utilisation HemoCue Hb 201*
2. Mode d'emploi (manuel d'utilisation) HemoCue Hb 201 DM
3. Dacie and Lewis, Practical Haematology, Ninth Edition, 2001
4. Reference and Selected Procedures for the Quantitative Determination of Hemoglobin in Blood; Approved Standard; CLSI Document H15-A
5. Interference Testing in Clinical Chemistry; Approved Guideline; CLSI Document EP7
6. Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline; CLSI Document EP5-A

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
1 = sangue venoso in EDTA, studio multicentrico
 2 = sangue capillare

Bibliografia

1. Manuale d'uso HemoCue Hb 201*
2. Istruzioni per l'uso di HemoCue Hb 201 DM (Manuale d'uso)
3. Dacie and Lewis, Practical Haematology, Ninth Edition, 2001
4. Reference and Selected Procedures for the Quantitative Determination of Hemoglobin in Blood; Approved Standard; CLSI Document H15-A
5. Interference Testing in Clinical Chemistry; Approved Guideline; CLSI Document EP7
6. Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline; CLSI Document EP5-A

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Literatur

1. HemoCue Hb 201+ Bedienungsanleitung
2. HemoCue Hb 201 DM Gebrauchsanweisung (Bedienungsanleitung)
3. Dacie and Lewis, Practical Haematology, Ninth Edition, 2001
4. Reference and Selected Procedures for the Quantitative Determination of Hemoglobin in Blood; Approved Standard; CLSI Document H15-A
5. Interference Testing in Clinical Chemistry; Approved Guideline; CLSI Document EP7
6. Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline; CLSI Document EP5-A

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+46 77 570 02 10
+46 77 570 02 12
info@hemocue.se
hemocue.com

151712 160127 GB

Expected Values
For general guidance to normal reference values applicable to most healthy adults and children, see Dacie and Lewis, Practical Haematology.
Due to a wide range of conditions (dietary, geographical, etc) which affect normal values, it is recommended that each laboratory establish its own normal range.

Specific Performance Characteristics

Repeatability and Within-laboratory Precision
Results summarized below was determined according to the CLSI Document EP5-A2.

Control level	N	x̄	SD	CV	SD	CV	Within-laboratory Precision	
							g/dL	%
1	400	7.70	0.057	0.74	0.100	1.30	0.71	
2	400	15.33	0.078	0.51	0.109	0.71		

Accuracy

Results of comparison studies between the HemoCue Hb 201 system and the International Council for Standardization in Haematology method (ICSH) are summarized below.

Study	N	Min	Max	Regression line	Correlation coefficient (r)
1	498	4.1	20.5	$Y = 1.009X - 0.008$	0.998
2	103	9.0	17.7	$Y = 0.971X + 0.584$	0.934

1 = venous EDTA blood, multicenter study
2 = capillary blood

Bibliography

- HemoCue Hb 201+ Operating Manual
- HemoCue Hb 201 DM Instructions for use (Operating Manual)
- Dacie and Lewis, Practical Haematology, Ninth Edition, 2001
- Reference and Selected Procedures for the Quantitative Determination of Hemoglobin in Blood; Approved Standard; CLSI Document H15-A
- Interference Testing in Clinical Chemistry; Approved Guideline; CLSI Document EP7
- Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline; CLSI Document EP5-A

Symbols used

	Caution		Batch code
	<i>In vitro</i> diagnostic medical device		Catalog number
	Do not reuse		Temperature limitation
	CE mark		Consult instructions for use
	Manufacturer		This product is covered by one or more patents
	Use by (Year Month Day)		Date of opening
Once opened – Expiry date		Open vial expiry date: Must not exceed the "Use by" date	

For markets outside the US and Canada only. For US and Canadian markets see last page.

HemoCue Hb 201 Microcuvettes, 4x25 & 4x50 Microcuvettes

The system and the HemoCue Hb 201 DM Analyzer (hereafter both types of analyzers are referred to as the HemoCue Hb 201 Analyzer).

Intended Use

Quantitative determination of hemoglobin in capillary, venous or arterial blood. The HemoCue Hb 201 systems are used for quantitative determination of hemoglobin in whole blood using a specially designed analyzer, the HemoCue Hb 201 Analyzer, and specially designed microcuvettes, the HemoCue Hb 201 Microcuvettes. HemoCue Hb 201 Microcuvettes are for *In Vitro* Diagnostic use only. The HemoCue Hb 201 Analyzer is only to be used with HemoCue Hb 201 Microcuvettes.

IVD Medical Device Directive

The HemoCue Hb 201 Microcuvettes comply with the IVD Medical Device Directive 98/79/EC and carry the CE mark.

Composition

The microcuvettes are made of polystyrene plastic. Reagents in µg/g microcuvette; <600 sodium deoxycholate, <300 sodium azide, <300 sodium nitrite, <350 nonreactive ingredients.

Precautions

Blood may be infectious, always handle with care and wear protective gloves. Consult local environmental authorities for proper disposal. For system information regarding precautions see Operating Manual.

Storage and Handling

Use HemoCue Hb 201 Microcuvettes prior to expiry date on each package. Microcuvettes in vial and individually packaged

Store at room temperature (15–30 °C, 59–86 °F) in a dry place. Once the seal of the vial is broken the microcuvettes in the vial are stable for three months.

Note: Always keep the vial properly closed.

Procedure

Operating temperature for HemoCue Hb 201+ System is 15–30 °C (59–86 °F) and operating temperature for HemoCue Hb 201 DM System is 18–30 °C (64–86 °F).

For proper use of the system, read Operating Manual regarding specimen collection, measurement/test procedure, maintenance, performance/theory and technical specifications.

Note: Make sure the microcuvette is completely filled before measurement

Quality Control

The HemoCue Hb 201 Analyzer has an internal electronic self-test, see Operating Manual. If a quality control test is required by local or other regulations, only use controls recommended by HemoCue. For further information regarding controls contact HemoCue AB.

Limitations

- a) Fill the microcuvette within 5 minutes after removal from the package.
- b) Start measurement no later than 10 minutes after filling the microcuvette.
- c) Confirm values above 23.5 g/dL (235 g/L, 14.6 mmol/L) with suitable laboratory method.
- d) Sulfhemoglobin is not measured with this method.
- e) Following substances have not been found to interfere: Acetaminophen (20 mg/dL), ascorbic acid (3 mg/dL), conjugated bilirubin (40 mg/dL), unconjugated bilirubin (20 mg/dL), creatinine (30 mg/dL), ibuprofen (40 mg/dL), leukocytes (600 x 10⁶ /L), lipemia (intralipid 4000 mg/L, triglycerides approximately 1000 mg/dL), salicylic acid (50 mg/dL), tetracycline (20 mg/dL), thrombocytes (2100 x 10⁹ /L), urea (500 mg/dL), uric acid (20 mg/dL). The highest concentration or percentages tested is referred to in brackets. Interference studies have been performed according to CLSI Document EP7.
- f) pH values between 6.3–9.0 do not interfere with the system.
- g) Mixing samples for an extended period can produce increased oxygen pressure and viscosity that may give falsely results.

151712 160127