



Afinion™ HbA1c External Quality Control Log

Afinion™ HbA1c Controls are recommended for routine quality control testing and should be analyzed:

- With each new shipment of Afinion™ test kits
- With each new lot of Afinion™ test kits
- At least every 30 days
- When training new operators in correct use of the Afinion™ HbA1c and the Alere Afinion™ AS100 Analyzer or the Afinion™ 2 Analyzer
- Anytime an unexpected test result is obtained
- If local, state and/or federal regulations require more frequent testing of control materials

Afinion™ Serial #: _____

AFINION™ HbA1c EXTERNAL CONTROLS ACCEPTABLE RANGES

	CONTROL C I	CONTROL C II
CONTROL LOT #/EXP		
OPEN VIAL EXPIRATION DATE		
% HbA1c CONTROL RANGE		

DATE	HbA1c TEST CARTRIDGE LOT#	HbA1c TEST CARTRIDGE EXPIRATION	CONTROL ID	% HbA1c RESULT	EXTERNAL CONTROL PASS/FAIL	TESTER'S INITIALS

REVIEWED BY: _____

DATE: _____

Abbott Diagnostics Technologies AS
 Kjelsåsveien 161
 P.O. Box 6863 Rodeløkka
 NO-0504 Oslo, Norway

Product support: Call 1.866.216.9505 | Abbott.com/poct