



PreView[®]

Strep A Test Kit




- Easy to Use
- FDA Approved and CLIA Waived
- Individually wrapped test strips
- Preliminary Results While Patient Waits
- External Control Included
- No Refrigeration Required
- Revenue Stream

PreView Diagnostics
— www.previewdiagnostics.com —


Test Procedure

1 Specimen Collection




2 Step 1

Add 4 drops Add 4 drops




3 Step 2



5X

1 MIN


4 Step 3




Marker Line

5 Step 4

5-10 MIN



6 Results



Positive Negative Invalid

Control line → C Test line → T

Order Information

Item Code	Description	Complexity
PRE-STA-25	25-Strip; Strep A Throat Swab Test	CLIA WAIVED