

DD-300 & DD-301 Non-display Doppler Systems



Newman  Medical

Employee Owned & Operated

800-267-5549

WWW.NEWMAN-MEDICAL.COM

You Deserve a Better Doppler!

Newman Medical offers a full line of Doppler systems with the performance and durability you expect to meet your application - obstetrical or vascular. The **DigiDop 300 and 301 Non-display Digital Doppler Systems** have unique features not found in other Doppler systems such as Digital Sound Clarity™. The DigiDop 301 adds the Smart Recharge System™. Each DigiDop is designed with the sensitivity you need while also designed to be tough and last through years of everyday office use. All our products are **Made in the USA in Colorado**. Each Newman Medical Doppler and probe must pass a series of comprehensive tests to assure quality and the durability of each product.

Probe choices include:

- ✦ 2MHz – obstetrical probe – later term
- ✦ 3MHz – obstetrical probe – general purpose
- ✦ 5MHz – vascular probe – deeper vessels
- ✦ 8MHz – vascular probe – general purpose
- ✦ 2MHz Waterproof – obstetrical water applications
- ✦ 3MHz Waterproof – obstetrical water applications
- ✦ Audio PPG – vascular PPG digit probe/sensor

Features:

- ✦ **Better Sound** – Digital Sound Clarity-DSC™ produces up to 6x signal of other Dopplers
- ✦ **Better Warranty** – 5 year honest warranty – all parts & labor. Includes rechargeable batteries on rechargeable units.
- ✦ **Better Support & Service** – Knowledgeable Staff & Guaranteed Service
- ✦ **Smart Recharge System™** (Model 301) A significant reduction in lifetime battery costs and also a reduction in waste
- ✦ 'Unbreakable' magnetic probe holder
- ✦ All probes are interchangeable
- ✦ User replaceable cord eliminates service downtime
- ✦ Powered by 3 AA batteries offer extended battery life of 900 mins
- ✦ Dimensions (h x w x l): 150 x 65 x 35 mm (6.0 x 2.5 x 1.25 in)
- ✦ Weight: 340 grams (12 ounces)
- ✦ Shipping Information: 2 pounds in 12"x9"x3" box