

Clarity Mono Cassettes



Product Code:
DTG-MONO (15/bx)

Contents:

- Mono Test Cassettes
- Developer Solution
- Sample Transfer Tubes
- Positive Control
- Negative Control
- Instructional Insert
- Procedure Card



- | | |
|---|---|
| <ul style="list-style-type: none">✓ Sensitivity >99%✓ Specificity 98.6%✓ 99% Overall Accuracy✓ Results as Early as 5 Minutes✓ No Age Restriction | <ul style="list-style-type: none">✓ 18 Month Dating✓ Room Temperature Storage✓ Test Can Run with Capillary Finger-Stick Whole Blood, Plasma, or Serum |
|---|---|



Distributed by STAT Technologies

For more information, contact
800-217-7828 or info@stat-technologies.com

www.stat-technologies.com

**We Support
Better Outcomes**

