

Clarity Strep A



Product Codes:
DTG-STP25 (25/bx)

Contents:

- Strep A Test Strips in a Vial
- Disposable Extraction Tubes
- Reagents A & B
- Positive Control
- Negative Control
- Workstation
- Procedure Card
- Package Insert

✓ Sensitivity: 95%

✓ Specificity: 98.4%

✓ 24 Month Dating Open Vial Stability

✓ Room Temperature Storage

✓ Controls Included in Kit

✓ Strips in Vial Format

✓ Made in the USA



 **Distributed by STAT Technologies**

For more information, contact
800-217-7828 or info@stat-technologies.com

www.stat-technologies.com

**We Support
Better Outcomes**