



URINE DRUG SCREENING DEVICE

iCUP™ Dx PRO 2

Detect opiates at low cutoff levels (300 ng/mL) and screen for buprenorphine

- Screen for 5-12 prescription and illicit drugs
- Simple to use, no need to tilt cup to activate test—read results at five minutes
- Self-contained cup ideal for sending preliminary positive specimens to the lab for confirmation
- Available with specimen validity test (SVT) and temperature strip—helps to monitor specimen tampering



For professional in vitro diagnostic use

iCup™ Dx Pro 2 Drugs of Abuse Cup

AVAILABLE DRUGS

Amphetamine (AMP)	Methadone (MTD)
Barbiturates (BAR)	Opiate (OPI300, OPI2000)
Buprenorphine (BUP)	Oxycodone (OXY)
Benzodiazepines (BZO)	Phencyclidine (PCP)
Cocaine (COC)	Tricyclic Antidepressants (TCA)
Methylenedioxy-methamphetamine (MDMA)	Marijuana (THC)
Methamphetamine (MET)	

SPECIMEN VALIDITY TESTS: Creatinine (CR), Glutaraldehyde (GL), Nitrites (NI), Oxidants (OX), Potential of Hydrogen (PH) and Specific Gravity (SG).

SIMPLE COLLECTION AND TESTING PROCEDURE

**iCup™ Dx
Device Components**

Cup lid

Peel off label

Results legend

Minimum fill line

Temperature strip

1 Collect specimen and secure the lid.

2 Peel label to reveal drug test and specimen validity.

3 Compare adulteration results to color chart. Interpret drug test results at 5 minutes.

PRODUCT PROCEDURE NOTE: Refer to product insert for complete instructions, limitations, precautions, and warnings. For *in vitro* diagnostic use at the point of care site. This device is a one-step, lateral-flow immunoassay device and provides only a preliminary qualitative test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography-mass spectrometry (GC-MS) or liquid chromatography-tandem mass spectrometry (LC-MS/MS) are the preferred confirmatory methods.

**TO LEARN MORE, CALL 800-340-4029
OR VISIT ABBOTT.COM/TOXICOLOGY**

© 2020 Abbott. All rights reserved. All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners. MKT50450 REV4 12/20

