




EXTEND YOUR MOBILITY

BRING POINT OF CARE TESTING CLOSER TO YOUR PATIENTS



AFINION™ 2 BATTERY KIT


Provides up to 10 hours (or approximately 60 tests) of continuous power for Afinion™ 2 analyzer, extending point-of-care testing for:

- mobile units
- home care visits
- drive-through testing sites
- nursing homes




AFINION™ 2 ANALYZER

Accurate results in minutes at every point-of-care location



AFINION™ 2 CARRY BAG

Can be worn on the back or shoulder for easy transport of Afinion™ 2 analyzer in support of on-the-go patient consultations



ORDERING INFORMATION

Product Name

Product Code

NEW AFINION™ 2 BATTERY KIT

1117292 (EMEA, EXCLUDING UK AND NORDICS)

NEW AFINION™ 2 BATTERY KIT

1117291 (UK)

AFINION™ CARRY BAG

B5592 (EMEA ONLY)

AFINION™ 2 ANALYZER

1116679

AFINION™ HBA1C

1116795

AFINION™ ACR

1116781

AFINION™ CRP

1116787

AFINION™ LIPID PANEL

1116801

AFINION™ 2 & CARRY BAG:

- Size of analyzer:
20.0 cm W x 18.6 cm H x 32.8 cm D
- Weight of analyzer: 3.40 kg
- Size of bag:
33.0 cm W x 19.0 cm H x 33.0 cm D
- Do not move the analyzer when a test cartridge is being processed
- Do not operate the analyzer in a moving vehicle
- Power down the analyzer before unplugging it

BATTERY KIT:

- Size of battery:
9.4 cm W x 6.8 cm H x 16.0 cm D
- Weight of battery: 1.0 kg
- Battery expected to last 10 hours or 60 tests (whichever comes first) before needing recharged
- Full charge expected to take up to 10 hours
- Use designated charger to charge battery
- Fits in carry bag



CONTACT SALES
FOR MORE
INFORMATION
(UK ONLY).



Watch this video to see how the Afinion™ 2 Carry Bag can be worn comfortably on your back or shoulder.