



TAP[®] Micro Select

Product Information Sheet

The TAP[®] Micro Select is a single use device designed to enable the sampling of whole capillary blood from the upper arm.

- FDA 510(k) Cleared
- CE Marked
- HALO[™] Technology



VIRTUALLY PAINLESS BLOOD SAMPLING

Our TAP[®] devices produce high-quality, whole blood samples at consistently high volumes of up to **900µL** through push-button activated, virtually painless capillary blood sampling.

Collect Up To
900µL

Whole Capillary Blood

Our HALO[™] Technology enables users to sample up to 900 µL of blood, accommodating a broad range of standard blood tests¹.

98%

Would Repeat Collection

When asked if they would repeat their sampling experience with the TAP[®] Micro Select device, 98% of patients responded "yes"¹.

¹ YourBio Health Data on File.



YourBioHealth.com

800-217-7828
info@stat-technologies.com

DEVICE SPECIFICATIONS

Packaged Dimensions: 122 x 117 x 58mm


Packaged Weight: ~30 grams

Sterilization Method: Gamma radiation

Storage: 18°C - 28°C (65°F - 82°F)

Shelf Life: 24 months

Manufactured in: Medford, MA, USA



HALO™ TECHNOLOGY

- Bladeless microneedle array
- Microneedles thinner than an eyelash
- Virtually painless compared to fingerstick, venipuncture, and competitive products

[STAT-Technologies.com](https://www.stat-technologies.com)

800-217-7828
info@stat-technologies.com

Handling

- Do not use if the package has been punctured.
- Do not use after expiration date.
- Single use device. Do not reuse.
- Do not resterilize.

Disposal

Dispose of device in accordance with local household sharps disposal laws and regulations. Packaging may be disposed of per local recycling guidelines or with common household waste.

Not made with
natural latex rubber

[YOURBIOHEALTH.COM](https://www.yourbiohealth.com)

© 2023 YourBio Health, Inc. All rights reserved.
TAP is a registered trademark of YourBio Health, Inc. The YourBio logo is a trademark of YourBio Health, Inc.

ISO
13485

 yourbio