



STAT Technologies rental programs are designed to provide additional testing equipment according to your needs without requiring the use of your capital budget.

For screening events that require additional equipment, STAT Technologies will provide the CardioChek Analyzers for a 15 to 90-day period at no charge with a qualifying test purchase. In addition, we offer a free 1-year Reagent Rental program designed for customers who need CardioChek analyzers for an extended period of time. In conjunction with PTS, we offer a no charge CardioChek System Placement Program for a minimum of 2 years with annual renewals available.

	CardioChek Loaner Program	CardioChek Rental Program	CardioChek System Placement Program
Terms:	<i>15 to 90 days*</i>	<i>1 Year w/ Annual Renewals</i>	<i>2 Years w/ Annual Renewals</i>
Test Type:	<i>All</i>	<i>All</i>	<i>Lipid Profile, Lipid+Glucose, or TC-HDL-Glucose</i>
Equipment Provided:	<i>CardioChek Analyzer</i>	<i>CardioChek Analyzer Only</i>	<i>CardioChek System</i>
Equipment Condition:	<i>Used</i>	<i>Refurbished</i>	<i>New</i>
Deposit Required:	<i>None</i>	<i>\$250.00 (Refundable)</i>	<i>None</i>
Cassette Commitment:	<i>\$300.00 (per analyzer per 15-day period)</i>	<i>\$500.00 (per analyzer per month)</i>	<i>4 bx Lipid Profile or Lipid+Glu 6 bx TC-HDL-Glucose</i>
Order Type:	<i>Customer Initiated</i>	<i>Customer Initiated</i>	<i>Monthly Auto-Ship</i>
Exclusive Purchase:	<i>Yes</i>	<i>Yes</i>	<i>Yes</i>
Cassette Deficiency Fees:	<i>\$75.00 (per analyzer per 15-day)</i>	<i>\$150.00 (per analyzer per month)</i>	<i>None</i>

*Terms include shipping to and from destination



CardioChek System – CLIA WAIVED

The CardioChek System empowers people to lead longer, healthier lives by providing easy to use, efficient, and economical diagnostic tools and information for immediate risk assessment and therapeutic monitoring of heart disease.

- Lab-Accurate—Certified by CDC’s and CRMLN programs and meets NCEP guidelines. FCC and ISO certified wireless communications
- Rapid Results—Results in 2 minutes using a simple 3-step process
- Reimbursable—Medicare and Private Insurance

1-800-217-7828

8760 Jefferson Hwy
Maple Grove, MN 55369



The PTS Panels® Lipid Panel+eGLU Smart Bundle 2-Pack (2729)
Includes 15 lipid panel test strips and 15 eGLU test strips



Three -Analyte



Lipid Panel Test Strip (1710)
Measures total cholesterol, HDL cholesterol and triglycerides
Calculates LDL cholesterol, total cholesterol/HDL ratio, non-HDL cholesterol



CHOL+HDL+GLU Panel Test Strip (2412)
Measures total cholesterol, HDL cholesterol and glucose
Calculates total cholesterol/HDL cholesterol ratio, non-HDL cholesterol

Use PTS Collect™ 40 µL capillary tubes (2866) with three-analyte test strips.

Dual -Analyte



CHOL+GLU Panel Test Strip (1765)
Measures total cholesterol and glucose



CHOL+HDL Panel Test Strip (1821)
Measures total cholesterol and HDL cholesterol
Calculates total cholesterol/HDL cholesterol ratio

Use PTS Collect 30 µL capillary tubes (2865) with dual-analyte test strips.

Single -Analyte



TOTAL CHOLESTEROL Test Strip (1711)
Measures total cholesterol



HDL CHOLESTEROL Test Strip (1714)
Measures HDL cholesterol



NO LONGER AVAILABLE - Triglyceride Tests



GLUCOSE Test Strip (1713)
Measures glucose



eGLU Test Strip (2713)
Measures electrochemical glucose

*Use PTS Collect 15 µL capillary tubes (2863) with single-analyte test strips.
No capillary tube required with eGLU test strips.*



Muti-Chemistry Controls (0721)



HDL Cholesterol Controls (0722)