

Clarity Drug Testing Cups



- ✓ Built-In Temperature Strip
- ✓ Results in 5 Minutes
- ✓ Screw Top Lid for a No-Leak Cup
- ✓ Multiple Drug Combinations Available
- ✓ 510K Cleared and CLIA Waived



DRUG OF ABUSE	ABBREVIATION	CUTOFF
Amphetamine	AMP	1000 ng/ml
Barbiturates	BAR	300 ng/ml
Buprenorphine/ Norbuprenorphine	BUP/NBUP	10 ng/ml
Benzodiazepine	BZO/BZD	300 ng/ml
Cocaine	COC	300 ng/ml
Methamphetamine	mAMP/MET	1000 ng/ml
MDMA or Ecstasy	MDMA	500 ng/ml
Morphine	MOP	300 ng/ml
Methadone	MTD	300 ng/ml
Opiate	OPI	2000 ng/ml
Oxycodone	OXY	100 ng/ml
Phencyclidine	PCP	25 ng/ml
Propoxyphene	PPX	300 ng/ml
Tricyclics	TCA	1000 ng/ml
Marijuana/Hashish	THC	50 ng/ml

Distributed by STAT Technologies

For more information, contact
800-217-7828 or info@stat-technologies.com

www.stat-technologies.com

**We Support
Better Outcomes**