



## TOXICOLOGY

FENTANYL  
1 ng/mL

OXYCODONE  
100 ng/mL

AMPHETAMINE  
500 ng/mL

MARIJUANA  
50 ng/mL

# COMPREHENSIVE CLIA-WAIVED<sup>1</sup> DRUG SCREENING

NOW WITH FENTANYL

USING THE iSCREEN™ URINE TEST DX  
DRUG SCREEN TOX CUP FROM THE  
GLOBAL LEADER IN TOXICOLOGY TESTING

With drug use at an all-time high and the U.S. opioid crisis declared a national emergency,<sup>2</sup> it is now more critical than ever for drug testing programs to have comprehensive screening menus that include fentanyl.

Be sure to include Abbott's new CLIA-waived<sup>1</sup> iScreen™ Urine Test DX Drug Screen Tox Cup in your testing program — now with fentanyl at detection levels as low as 1 ng/mL. This cup also includes 13 other high prevalence substances to help keep your testing program and participants accountable.



## WHY THIS MATTERS: TESTING NEEDS TO CONTINUE TO EVOLVE

In 2022, fentanyl was a contributing factor to 90% of opioid-related deaths. Fentanyl is 50x more potent than heroin and 100x more potent than morphine and can unknowingly be laced with other drugs, making them cheaper, more powerful, more addictive and more dangerous. Many people do not even know they are taking fentanyl until it's too late. Now is the time to include fentanyl in your testing program to help identify and stop its harmful use.

## OUR SOLUTION: A NEW CLIA-WAIVED<sup>1</sup> CUP WITH FENTANYL

A simple-to-use, rapid-result drug screen focused on testing for prevalent drugs at low cutoffs, including:

CLIA-WAIVED <sup>1</sup> 14-PANEL iSCREEN™ URINE TEST DX DRUG SCREEN TOX CUP – URINE		
Drug	Abbreviation	Cutoff (ng/mL)
Amphetamine	AMP	500
Barbiturates	BAR	300
Benzodiazepines	BZO	300
Buprenorphine	BUP	10
Cocaine	COC	150
Ecstasy	MDMA	500
Fentanyl	FTY	1
Marijuana	THC	50
Methadone	MTD	300
Methamphetamine	MET	500
Opiates/Morphine	MOP	300
Oxycodone	OXY	100
Phencyclidine	PCP	25
Tricyclic Antidepressant	TCA	1,000

Also Includes Specimen Validity Testing of Creatinine, pH and Oxidants

## WHY CHOOSE THE iSCREEN™ URINE TEST DX DRUG SCREEN TOX CUP?

Support your drug testing program with the iScreen™ Urine Test DX Drug Screen Tox Cup, and take the next step in promoting healthier and safer communities.

### FAST AND ACCURATE

Offer accurate results in as little as 5 minutes, with results remaining stable for up to 10 minutes

### CLIA-WAIVED<sup>1</sup> AND EASY TO USE

Implement a simple, one-step solution with a low risk for errors (FDA 510(k) cleared and CLIA-waived<sup>1</sup>)

### BROAD MENU

Test for 14 high prevalence drugs, including fentanyl, while also testing for sample adulteration

Contact us today to learn more about how the iScreen™ Urine Test DX Drug Screen Tox Cup can support your initiatives for healthier and safer communities.

1. Indications and clearance: 510(k) cleared and CLIA-waived for over the counter (OTC) use.
2. Renewal of Determination that a Public Health Emergency Exists Nationwide | Opioid-September 26, 2023. <https://aspr.hhs.gov/legal/PHE/Pages/Opioid-26Sept2023.aspx#:~:text=%C2%A7%20247d%2C%20do%20hereby%20renew>

© 2024 Abbott. All rights reserved. All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners. Any photos displayed are for illustrative purposes only. COL-26651 v1 11/24

DISTRIBUTED BY STAT TECHNOLOGIES



800-217-7828

INFO@STAT-TECHNOLOGIES.COM

WWW.STAT-TECHNOLOGIES.COM