

REDDROP ONE SPEC SHEET

Spec Sheet

Regulatory Status	FDA 510k Cleared
Product Type	Lancet with permanent microneedle retraction
Patient Segment	No age restrictions on clearance
Disposable	Yes
HSA/FSA Eligible	Yes
Microneedle Length	1.75mm
Minimum Order Quantity	80
Sterilized	Yes
Sterilization Method	Gamma Radiation
Storage	18°C to 28°C (64.4°F to 82.4°F)
Shelf Life	12 Months
BD Microtainer Compatible	Yes
Packaged Weight	42 Grams



- ✓ No needles
- ✓ No Tears
- ✓ No missed veins
- ✓ 99% Satisfaction Rate