


accutest[®]

ORAL FLUID COTININE TEST

CLIA COMPLEXITY:
WAIVED

CAT#DS443

RESULTS INTERPRETATION



Saliva Cotinine Test

The Accutest Saliva Cotinine Test is a rapid, screening test for the simultaneous, qualitative detection of Cotinine in human oral fluid. The test detects Cotinine in oral fluids at a cutoff of 30ng/mL.


This product can be used by CLIA Waived facilities.

FEATURES & BENEFITS

- OTC - Over-the-Counter ✓
- CLIA Categorized: "Waived" ✓
- Rapid, On-Site Determination ✓
- Simple Collection ✓
- Results in 10 minutes ✓



Simple Collection Procedure



- 1** Collect oral fluid.
- 2** Saturation indicator on collection swab turns red to confirm sufficient volume is collected.
- 3** Results in 10 min after collection swab is screwed into test device.

JANT PHARMACAL CORPORATION
LABORATORY & POINT OF CARE DIAGNOSTICS

accutest[®] is a registered trademark of Jant Pharmacal Corporation.

Jant Pharmacal Corporation
16530 Ventura Blvd.,
Suite 512 Encino, CA 91436
800.676.5565 Fax: 818.986.0235
info@JantDX.com JantDX.com

CONTACT US FOR MORE INFO
800.676.5565