

# DIGIDOP PROBES

FOR DIGIDOP CLASSIC AND DIGIDOP II DOPPLER SYSTEMS



American-Made Quality



Includes labor, parts, and accidental breakage

## FETAL



### D2

2MHz fetal probe - optimized for later term pregnancies or larger patients.



### D2W

2MHz waterproof fetal probe - designed for water applications and labor deliveries.



### D3

3MHz general purpose obstetrical probe - optimized for early fetal heartbeat detection and use throughout pregnancy.



### D3W

3MHz waterproof fetal probe - designed for water applications and labor deliveries.

## VASCULAR



### D5

5MHz probe - optimized for deeper vessels in vascular applications.



### D5-SD

5MHz SteadyDop probe - wide sensors for easy contact with deeper arteries.



### D8

8MHz general purpose vascular probe - optimized for peripheral vessels.



### D8-SD

8MHz SteadyDop probe - wide sensors for easy contact with peripheral arteries.



### DPPG

PPG probe for vascular applications - use when Doppler probe is insufficient (incompressible arteries).