



Value+TM One-step hCG Tests



For the Early Detection of Pregnancy

Accurate Results

- > 99% Sensitivity
- 25 mIU/mL in urine or serum
- Built-in procedural control

Easy to Use

- One-step procedure

Rapid Results

- 5 minutes or less

**CLIA
Waived**

Convenient

- Room temperature storage
- 24 months shelf life (from date of manufacture)

Made in the USA 

Accutest[®] Value+TM Pregnancy Test Device

hCG URINE CASSETTE

- 25 tests per kit
- Catalog No. PF853
- CLIA Waived
- CPT Code: 81025

hCG SERUM/URINE CASSETTE


- 25 tests per kit
- Catalog No. PF864
- CLIA Waived for Urine
- CPT Code: 84703



Accustrip[®] Value+TM hCG Test Strip

hCG URINE STRIP

- 25 tests per kit
- Catalog No. PF851
- CLIA Waived
- CPT Code: 81025



accutest[®]

JANT PHARMACAL CORPORATION

16530 Ventura Boulevard, Suite 512, Encino, California 91436

800.676.5565 818.986.8530 Fax 818.986.0235

www.accutest.net info@accutest.net Accutest and Accustrip are registered trademarks of Jant Pharmacal Corporation M018AC-1

