


LH Ovulation Rapid Test Strip



Intended Use

Primarily used to detect the level of luteinizing hormone (LH). Convenient for calculating ovulation period and planning for pregnancy


Specifications

- ✓ **Type of Use:** OTC
- ✓ **Sample Type:** Urine (one-time use urine cup provided)
- ✓ **Test Format:** Test Strip
- ✓ **Cut-off Levels:** 25 mIU/mL
- ✓ **Results Time:** 3-5 minutes
- ✓ **Packaging:** 50 strips (customizable)
- ✓ **Qualification:** FDA 510(k) Exempt, CE
- ✓ **Wider Strip Design:** 5mm
- ✓ **Over 99% accuracy**

AI-Interpreted Ovulation Test Result



<https://pregnancynow.ai>



Problem Addressed

Interpretation Challenges: LH test results require users to observe and interpret them on their own. Sometimes, the color changes may not be very apparent, which can lead to misinterpretation. Based on the AI intelligence program developed for ovulation test strips use with the test strips.

Core Functions AI image recognition: Quantify LH levels, accurately tracking and predicting ovulation periods, helping users maximize conception opportunities and avoid misjudgment. In addition, with user's menstrual period information, send reminders for testing times, and ovulation predictions.

Simplify Operation: Simplify the ovulation detection process and reduce cumbersome steps through features such as test strip scanning.