

ZOLL AED 3 Smart. Ready. Trusted



Enhanced Real CPR Help® — our proven real-time CPR feedback technology now features a full-color display with vivid rescue images, a CPR cycle timer, and a large color bar gauge that shows CPR compression depth.

Innovative Design and Enhanced Features

The ZOLL AED 3® defibrillator builds on our reputation for providing rescuers with unrivaled rescue support and a low total cost of ownership. The ZOLL AED 3 embodies all of the same great qualities of the AED Plus® defibrillator, while offering new, advanced features, including:

- **Real CPR Help** — Provides integrated, real-time feedback on CPR compression rate and depth.
- **Integrated pediatric rescue** — Universal CPR Uni-padz® electrodes can be used to treat both adult and pediatric victims of sudden cardiac arrest (SCA) by simply activating child mode.
- **WiFi connectivity** — Automatically reports AED readiness.
- **Long-life consumables** — Five-year battery and pads reduce AED device maintenance requirements and the total cost of ownership of your ZOLL AED 3.

ZOLL AED 3 Package:

Package includes ZOLL AED 3, CPR Uni-padz (for adult and pediatric use), ZOLL AED 3 battery pack, and manuals. Free on-site in-service with every purchase.

Semi-automatic ZOLL AED 3 Package Part# 8511-001101-01

Automatic ZOLL AED 3 Package Part# 8511-001102-01



350,000

Out-of-hospital cardiac arrests in 2016¹



10%

Decrease in survival for every minute of delayed defibrillation²



5%

Chance of survival if a victim collapses and nothing is done until emergency services arrive³

¹American Heart Association website. http://cpr.heart.org/AHA/ECC/CPRandECC/General/UCM_477263_Cardiac-Arrest-Statistics.jsp. Accessed 21 July 2020.

²Weisfeldt ML, et al. *J Am Coll Cardiol*. 2010;55(16):1713–20.

³*Circulation*. 2010;122:S720–S728, <https://doi.org/10.1161/CIRCULATIONAHA.110.970970>.