




SAFE life™

T-Dip Multi-Drug Urine Test Panel

Rapid Point-of-Care Toxicology Testing

- CLIA Waived, 510(K) Cleared & Forensic Use Only
- Cost Effective & Customizable
- Test up to 16 Drugs per Card



CLIA-waived Drugs & Cutoff Levels

AMP AMP AMP	Amphetamine (AMP)	1000, 500 ng/mL	TCA TCA TCA	Tricyclic Antidepressants (TCA)	1000 ng/mL
BAR BAR BAR	Barbiturates (BAR)	300 ng/mL	MET MET MET	Methamphetamine (MET/mAMP)	1000, 500 ng/mL
BZO BZO BZO	Benzodiazepines (BZO)	300 ng/mL	MDMA MDMA MDMA	Methylenedioxymethamphetamine (MDMA)	500 ng/mL
BUP BUP BUP	Buprenorphine (BUP)	10 ng/mL	MOP MOP MOP	Morphine (MOP/OPI300)	300 ng/mL
COC COC COC	Cocaine (COC)	300, 150 ng/mL	OPI OPI OPI	Opiate (OPI)	2000 ng/mL
THC THC THC	Cannabinoids (THC)	50 ng/mL	OXY OXY OXY	Oxycodone (OXY)	100 ng/mL
EDDP EDDP EDDP	2-ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine (EDDP)	300 ng/mL	PCP PCP PCP	Phencyclidine (PCP)	25 ng/mL
MTD MTD MTD	Methadone (MTD)	300 ng/mL	PPX PPX PPX	Propoxyphene (PPX)	300 ng/mL

Forensic Use Panels

6-MAM 6-MAM 6-MAM	Monoacetylmorphine (6-MAM)	10 ng/mL	K2 K2 K2	Synthetic Cannabinoid (K2)	50 ng/mL
AMP AMP AMP	Amphetamine (AMP)	300, 500 ng/mL	KET KET KET	Ketamine (KET)	300, 1000 ng/mL
BUP BUP BUP	Buprenorphine (BUP)	5 ng/mL	KRA KRA KRA	Kratom (KRA)	100, 300 ng/mL
BZO BZO BZO	Benzodiazepines (BZO)	100, 200 ng/mL	LSD LSD LSD	Lysergic Acid Diethylamide (LSD)	20 ng/mL
COC COC COC	Cocaine (COC)	100, 150 ng/mL	MET MET MET	Methamphetamine (MET/mAMP)	300, 500 ng/mL
COT COT COT	Cotinine (COT)	200 ng/mL	MOP MOP MOP	Morphine (MOP or OPI)	100 ng/mL
EDDP EDDP EDDP	EDDP	100 ng/mL	MTD MTD MTD	Methadone (MTD)	200 ng/mL
EtG EtG EtG	Ethyl Glucuronide (EtG)	500 ng/mL	MQL MQL MQL	Methaqualone (MQL)	300 ng/mL
FTY FTY FTY	Fentanyl (FTY)	20 ng/mL	PGB PGB PGB	Pregabalin (PGB)	500 ng/mL
GAB GAB GAB	Gabapentin (GAB)	2000 ng/mL	THC THC THC	Cannabinoids (THC)	15, 25, 40 ng/mL
HMO HMO HMO	Hydromorphone (HMO)	300 ng/mL	TRA TRA TRA	Tramadol (TRA)	100, 200, 1000 ng/mL

Generic stock adulteration testings include Creatinine (CR), Specific Gravity (S.G.) and pH. Other adulterates like Glutaraldehyde (GL), Oxidant (OX) & Nitrite (NI) are available through special order.



SAFE life™

T-Dip Single Drug Urine Test Card

※ Kit Size: 25 / Case Size: 20

Item Code	Description	Complexity
WDBU-114	1-Panel Card; BUP Buprenorphine	CLIA WAIVED
WDBZ-114	1-Panel Card; BZO Benzodiazepines	CLIA WAIVED
WDCO-114	1-Panel Card; COC Cocaine	CLIA WAIVED
WDMA-114	1-Panel Card; mAMP/MET Methamphetamine	CLIA WAIVED
WDMT-114	1-Panel Card; MTD Methadone	CLIA WAIVED
WDOP-114	1-Panel Card; OPI Opiate 2000 ng/mL	CLIA WAIVED
WDOX-114	1-Panel Card; OXY Oxycodone	CLIA WAIVED
WDPC-114	1-Panel Card; PCP Phencyclidine	CLIA WAIVED
WDTH-114	1-Panel Card; THC Marijuana/Cannabinoids	CLIA WAIVED
W6AM-114	1-Panel Card; 6-MAM Monoacetylmorphine	Forensic Use Only
WCOT-114	1-Panel Card; COT Cotinine	Forensic Use Only
WDK2-114	1-Panel Card; K2 Synthetic Cannabinoid	Forensic Use Only
** WDTH-15	1-Panel Card; THC Marijuana 15 ng/mL	Forensic Use Only
WDTR-114	1-Panel Card; TRA Tramadol	Forensic Use Only
WETG-114	1-Panel Card; EtG Ethyl Glucuronide	Forensic Use Only
WETG-300	1-Panel Card; EtG Ethyl Glucuronide 300 ng/mL	Forensic Use Only
WFTY-114	1-Panel Card; FTY Fentanyl	Forensic Use Only
** WGAB-114	1-Panel Card; GAB Gabapentin	Forensic Use Only
** WHMO-114	1-Panel Card; HMO Hydromorphone	Forensic Use Only
WKET-114	1-Panel Card; KET Ketamine 1000 ng/mL	Forensic Use Only
WKRA-114	1-Panel Card; KRA Kratom	Forensic Use Only
WMQL-114	1-Panel Card; MQL Methaqualone	Forensic Use Only
** WPGB-114	1-Panel Card; PGB Pregabalin	Forensic Use Only

** Ordered upon request

T-Dip Multi-Drug Urine Test Panel

※ Kit Size: 25 / Case Size: 20

Item Code	Description	Complexity
WDOA-234	3-Panel Card; COC,mAMP,THC	CLIA WAIVED
WDOA-154	5-Panel Card; AMP,COC,OPI,PCP,THC	CLIA WAIVED
WDOA-254	5-Panel Card; AMP,COC,mAMP,OPI,THC	CLIA WAIVED
WDOA-654	5-Panel Card; BZO,COC,mAMP,OPI,THC	CLIA WAIVED
WDOA-754	5-Panel Card; AMP,BZO,COC,OPI,THC	CLIA WAIVED
WDOA-264	6-Panel Card; AMP,BZO,COC,mAMP,OPI,THC	CLIA WAIVED
WDOA-564	6-Panel Card; BZO,COC,mAMP,OPI,OXY,THC	CLIA WAIVED
WDOA-865NT	6-Panel Card without THC; BUP,BZO,COC,MOP,MTD,OXY	CLIA WAIVED
WDOA-274	7-Panel Card; AMP,BZO,COC,mAMP,OPI,OXY,THC	CLIA WAIVED
WDOA-295	9-Panel Card; AMP,BUP,BZO,COC,MOP,MTD,OXY,PPX,THC	CLIA WAIVED
WDOA-994NT	9-Panel Card without THC; AMP,BAR,BUP,BZO,COC,MET,MTD,OPI,OXY	CLIA WAIVED
WDOA-1104	10-Panel Card; AMP,BAR,BZO,COC,mAMP,MTD,OPI,PCP,TCA,THC	CLIA WAIVED
WDOA-3104	10-Panel Card; AMP,BAR,BZO,COC,mAMP,MDMA,MTD,OPI,PCP,THC	CLIA WAIVED
WDOA-4104	10-Panel Card; AMP,BAR,BZO,COC,mAMP,MTD,OPI,OXY,PCP,THC	CLIA WAIVED
WDOA-7104	10-Panel Card; BAR,BUP,BZO,COC,mAMP,MTD,OPI,OXY,TCA,THC	CLIA WAIVED
WDOA-9104	10-Panel Card; AMP,BAR,BUP,BZO,COC,mAMP,MTD,OPI,OXY,THC	CLIA WAIVED
WDOA-1124	12-Panel Card; AMP,BAR,BZO,COC,mAMP,MDMA,MTD,OPI,OXY,PCP,PPX,THC	CLIA WAIVED
WDOA-3124	12-Panel Card; AMP,BAR,BUP,BZO,COC,MDMA,MTD,TCA,OPI,OXY,PCP,THC	CLIA WAIVED
WDOA-6124	12-Panel Card; AMP,BAR,BUP,BZO,COC,mAMP,MDMA,MTD,OPI,OXY,PCP,THC	CLIA WAIVED
WDOA-6125	12-Panel Card; AMP,BAR,BUP,BZO,COC,mAMP,MDMA,MOP,MTD,OXY,PCP,THC	CLIA WAIVED
WDOA-7125	12-Panel Card; AMP,BAR,BZO,COC,mAMP,MDMA,MOP,MTD,OXY,PCP,TCA,THC	CLIA WAIVED
WDOA-8145EF	14-Panel Card; AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, THC, ETG, FTY	Forensic Use Only

Rel.: 2024/05/10